

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

Elektrodendochte Standardgröße	Kat.Nr./Cat.No. :	42988
Electrode wicks standard size, Filter cardboard	Lot/Contr.No. :	100344

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Größe size		120 x 6 x 1mm	entspricht corresponds
Anwendung application		für Isoelektrische Fokussierung for Isoelectric focussing	
Mindesthaltbarkeit minimum shelf life			04/ 2015
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 30.06.2014

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.