

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

| | |
|---|--------------------------------|
| Elektrodenochte Standardgröße | Kat.Nr./Cat.No. : 42988 |
| Electrode wicks standard size, Filter cardboard | Lot/Contr.No. : 210498 |

| Parameter parameter | Methode method | Spezifikation specification | Ergebnis result |
|---|--------------------------|--|---------------------------|
| Größe size | | 120 x 6 x 1mm | entspricht corresponds |
| Anwendung application | | für Isoelektrische Fokussierung for Isoelectric focussing | |
| Mindesthaltbarkeit minimum shelf life | | | 27.11.2025 |
| Lagerung (°C) storage | | | +15 bis/ to +30 |

SERVA Electrophoresis GmbH

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 21.04.2021

Christian Monsler

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.