

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

GEL-FIX® for PAG 265 x 125 mm	Kat.Nr./Cat.No. : 42993 Lot/Contr.No. : P160303
--	--

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Beschreibung description	klare, farblose Haftfolien aus Polyester, auf beiden Seiten aktiviert mit einer Acrylatschicht clear, colorless polyester sheets, activated on both sides to bind polyacrylamide gels		
Größe (mm) size		265 x 125	entspricht corresponds
Dicke (mm) thickness		0.18	entspricht corresponds
Anwendung application	zum Herstellen von Polyacrylamidgelen for producing polyacrylamide gels		
Mindesthaltbarkeit minimum shelf life			13.07.2019
Lagerung (°C) storage			+15 bis/ to +30

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 12.09.2019

Christian Monsler

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.