

ANALYSENZERTIFIKAT  
CERTIFICATE OF ANALYSIS



<b>GEL-FIX® for PAG</b> <b>265 x 125 mm</b>	<b>Kat.Nr./Cat.No. : 42993</b> <b>Lot/Contr.No. : P230440</b>
--	--

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Beschreibung</b> description	klare, farblose Haftfolien aus Polyester, auf beiden Seiten aktiviert mit einer Acrylatschicht clear, colorless polyester sheets, activated on both sides to bind polyacrylamide gels		
<b>Größe (mm)</b> size		265 x 125	entspricht corresponds
<b>Dicke (mm)</b> thickness		0.18	entspricht corresponds
<b>Anwendung</b> application	zum Herstellen von Polyacrylamidgelen for producing polyacrylamide gels		
<b>Mindesthaltbarkeit</b> minimum shelf life			13.12.2026
<b>Lagerung (°C)</b> storage			+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 02.02.2024**

Dr. Judith Koch

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.