

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

GEL-FIX[®] for Covers	Kat.Nr./Cat.No. :	42995
	Lot/Contr.No. :	P110053

Parameter parameter	Spezifikation specification	Ergebnis result
Beschreibung description	klare, farblose Folien aus Polyester; auf beiden Seiten mit einer hydrophilen Schicht überzogen clear, colorless polyester sheets; hydrophilic layer to prevent from adhesion	
Größe (mm) size	280 x 125	entspricht corresponds
Dicke (mm) thickness	0.075	entspricht corresponds
Anwendung application	zum Abdecken von Gelen for covering gel surfaces	
Mindesthaltbarkeit minimum shelf life		02/ 2016
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 28.08.2012

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.