

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

<b>SERVAGel(TM) IEF Starter Kit</b>	<b>Kat.Nr./Cat.No. : 43205</b>
	<b>Lot/Contr.No. : P240287</b>

<b>Parameter</b> parameter	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Anwendung</b> application	Für die isoelektrische Fokussierung von Proteinen For the isoelectric focusing of proteins	entspricht corresponds
<b>Inhalt</b> content	SERVAGel(TM) IEF 3 – 10 SERVA IEF Anode Buffer, powder SERVA IEF Cathode Buffer, powder IEF sample buffer (2x) IEF Marker 3-10, Liquid Mix SERVA Violet 17	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life		05.06.2025
<b>Lagerung (°C)</b> storage		+2 bis/ to +8

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 20.01.2025**

Dr. Judith Koch

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.