

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

Buffer Kit for 2D HPE(TM) Gels	Kat.Nr./Cat.No.:	43312
	Lot/Contr.No.:	240946

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Beschreibung description	Kit enthält Anoden- und Kathodenpuffer, Äquilibrierungspuffer, Kühlkontaktflüssigkeit und Elektrodendochte. Kit comprising anode, cathode and equilibration buffer, cooling contact fluid and electrode wicks.		
Mindesthaltbarkeit minimum shelf life			21.03.2027
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 05.11.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.