

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

2D Gel DALTsix NF 12.5 % Kit	Kat.Nr./Cat.No.:	43313
	Lot/Contr.No.:	151142-151494

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Beschreibung description	Kit enthält 6 an Kunststoffolie gebundene Gele, einschließlich Anoden- und Kathodenpuffer, IPG Equilibrierungspuffer und Agarose. Format: 260 x 200 x 1 mm Kit comprising 12 plastic-backed gels, including anodal and cathodal buffer, IPG equilibration buffer and agarose. Size: 260 x 200 x 1 mm		
Anwendungstest application test	SDS-Elektrophorese SDS-Electrophoresis	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			06.05.2017
Lagerung (°C) storage	Buffer Kit 2D Gel DALTsix		+15 bis/ to +30
Lagerung (°C) storage	2D Gel Dalt NF 12.5%		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 18.12.2015

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.