

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Buffer Kit 2D Gel DALTsix	Kat.Nr./Cat.No.:	43325
	Lot/Contr.No.:	170298

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Beschreibung description	Kit enthält Anoden- und Kathodenpuffer, IPG Equilibrierungspuffer, Doppelklebeband und Agarose. Kit comprising anodal and cathodal buffer, IPG equilibration buffer, double adhesive tapeand and agarose.		
Mindesthaltbarkeit minimum shelf life			23.09.2018
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 28.03.2017

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.