

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

| | | |
|----------------------------------|-------------------------|---------------|
| Buffer Kit 2D Gel DALTsix | Kat.Nr./Cat.No.: | 43325 |
| | Lot/Contr.No.: | 171040 |

| Parameter parameter | Methode method | Spezifikation specification | Ergebnis result |
|---|--|---------------------------------------|---------------------------|
| Beschreibung description | Kit enthält Anoden- und Kathodenpuffer, IPG Equilibrierungspuffer, Doppelklebeband und Agarose. Kit comprising anodal and cathodal buffer, IPG equilibration buffer, double adhesive tapeand and agarose. | | |
| Mindesthaltbarkeit minimum shelf life | | | 31.08.2018 |
| Lagerung (°C) storage | | | +15 bis/ to +30 |

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 21.11.2017

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.