

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

Buffer Kit 2D Gel DALTsix	Kat.Nr./Cat.No.:	43325
	Lot/Contr.No.:	191187

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Beschreibung description	Kit enthält Anoden- und Kathodenpuffer, IPG Equilibrierungspuffer, Doppelklebeband und Agarose. Kit comprising anodal and cathodal buffer, IPG equilibration buffer, double adhesive tapeand and agarose.		
Mindesthaltbarkeit minimum shelf life			28.06.2020
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 04.12.2019

Christian Monsler

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.