

ANALYSENZERTIFIKAT
 CERTIFICATE OF ANALYSIS

IgG Sample Diluter IEF	Kat.Nr./Cat.No. :	43336
	Lot/Contr.No. :	240970

Parameter parameter	Spezifikation specification	Ergebnis result
Aussehen appearance	klare rote Lösung clear red solution	entspricht corresponds
Anwendung application	Verdünnungspuffer für CSF-Analyse in IEF-Gelen dilution buffer required for CSF analysis in IEF gels	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		21.10.2026
Lagerung (°C) storage		+2 bis/ to +8

SERVA Electrophoresis GmbH

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 11.11.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.