

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS



---

<b>EPO Doping IEF Kit 30S</b>	<b>Kat.Nr./Cat.No. :</b>	<b>43389</b>
	<b>Lot/Contr.No. :</b>	<b>250324</b>

---

<b>Parameter</b> parameter	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Anwendung</b> application	EPO IEF Analytik auf Polyacrylamidgelen EPO IEF analysis with polyacrylamide gels	entspricht corresponds
<b>Inhalt</b> content	Clean Gel EPO IEF 30S Buffer Kit For EPO Electrode Wicks	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life		20.05.2026
<b>Lagerung (°C)</b> storage		-15 bis/ to -25

SERVA Electrophoresis GmbH

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b> <b>Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.</b> We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date:** 02.04.2025

Dr. Judith Koch

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.