

ANALYSENZERTIFIKAT
 CERTIFICATE OF ANALYSIS

EPO Doping IEF Kit 30S	Kat.Nr./Cat.No. :	43389
	Lot/Contr.No. :	P140557-180517

Parameter parameter	Spezifikation specification	Ergebnis result
Anwendung application	EPO IEF Analytik auf Polyacrylamidgelen EPO IEF analysis with polyacrylamide gels	entspricht corresponds
Inhalt content	Clean Gel EPO IEF 30S Buffer Kit For EPO Electrode Wicks	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		09.12.2021
Lagerung (°C) storage		-15 bis/ to -25

SERVA Electrophoresis GmbH

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH Qualitätskontrolle/ Quality Control	Druckdatum/ printing date: 10.12.2020
---	--

Dr. Judith Koch Marina Konter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.