

ANALYSENZERTIFIKAT
 CERTIFICATE OF ANALYSIS

EPO Doping IEF Kit 30S	Kat.Nr./Cat.No. :	43389
	Lot/Contr.No. :	P210280-230295

Parameter parameter	Spezifikation specification	Ergebnis result
Anwendung application	EPO IEF Analytik auf Polyacrylamidgelen EPO IEF analysis with polyacrylamide gels	entspricht corresponds
Inhalt content	Clean Gel EPO IEF 30S Buffer Kit For EPO Electrode Wicks	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		05.01.2024
Lagerung (°C) storage		-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 13.03.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.