

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

CSF Analysis Kit for PhastSystem	Kat.Nr./Cat.No. :	43393
	Lot/Contr.No. :	140513

Parameter parameter	Spezifikation specification	Ergebnis result
Anwendung application	Auftrennung von CSF-Proben separation of CSF samples	entspricht corresponds
Inhalt content	CleanGel IEF for PhastSystem IEF Mix for CSF Analysis IEF sample application strips 8 slots	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		11/ 2015
Lagerung (°C) storage		-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 13.10.2014

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.