

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

SERVA HPE(TM) Silver Staining Kit	Kat.Nr./Cat.No. :	43395
	Lot/Contr.No. :	170155

Parameter parameter	Spezifikation specification	Ergebnis result
Inhalt content	Kit enthält Reagenzien zum Färben von 4 Large 2D HPE™ Gelen. Kit contains reagents to stain 4 large 2D HPE™ gels.	
Mindesthaltbarkeit minimum shelf life		05.08.2019
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control**Druckdatum/ printing date: 07.12.2017**

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.