

ANALYSENZERTIFIKAT

CERTIFICATE OF ANALYSIS

SERVA HPE™ Coomassie® Staining Kit	Kat.Nr./Cat.No. : 43396
	Lot/Contr.No. : P230452

Parameter parameter	Spezifikation specification	Ergebnis result
Anwendung application	Sensitive Kolloidalfärbung von 1D- und 2D-Gelen nach der Elektrophorese. Highly sensitive colloidal staining of 1D and 2D gels after electrophoresis.	entspricht corresponds
Inhalt content	1000 ml Solution A 25 ml Solution B	
Mindesthaltbarkeit minimum shelf life		15.12.2026
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
 Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 23.12.2025

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.