

ANALYSENZERTIFIKAT
 CERTIFICATE OF ANALYSIS

SERVA CSF Silver Staining Kit	Kat.Nr./Cat.No. : 43398
	Lot/Contr.No. : P240219

Parameter parameter	Spezifikation specification	Ergebnis result
Anwendung application	Hochsensitiver Silberfärbekit für das Färben von bis zu 5 foliengestützten, ultra-dünnen Horizontal-IEF-Gelen für die CSF-Analyse. High sensitive silver staining kit for staining of up to 5 film-backed, ultra-thin horizontal IEF gels for CSF analysis.	entspricht corresponds
Inhalt content	Fixing Solution I Ethanol denatured for Fixing Solution II and Wash Solution Fixing Solution III Solution A Solution B Solution C Solution D Solution E	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		31.03.2025
Lagerung (°C) storage		+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
 Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH Qualitätskontrolle/ Quality Control	Druckdatum/ printing date: 31.07.2024
---	--

Dr. Judith Koch	Daniela Lux-Helmstetter
-----------------	-------------------------

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.