

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS



---

<b>SERVA CSF Silver Staining Kit</b>	<b>Kat.Nr./Cat.No. :</b> 43398
	<b>Lot/Contr.No. :</b> P240309

---

<b>Parameter</b> parameter	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Anwendung</b> application	Hochsensitiver Silberfärbekit für das Färben von bis zu 5 foliengestützten, ultra-dünnen Horizontal-IEF-Gelen für die CSF-Analyse. High sensitive silver staining kit for staining of up to 5 film-backed, ultra-thin horizontal IEF gels for CSF analysis.	entspricht corresponds
<b>Inhalt</b> content	Fixing Solution I Ethanol denatured for Fixing Solution II and Wash Solution Fixing Solution III Solution A Solution B Solution C Solution D Solution E	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life		30.09.2025
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b> <b>Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.</b> We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
---

---

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date:** 14.10.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.