

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

IPG Chamber Cleaner	Kat.Nr./Cat.No. : 43399
IPG Chamber Cleaner	Lot/Contr.No. : P130057

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		gelbe, viskose Flüssigkeit yellow, viscous liquid	entspricht corresponds
Anwendung use		Reinigung cleaning procedures	
pH	1:10 verd./ dil. in H ₂ O	7.0 – 8.0	7.0
Leitfähigkeit (mS/ cm) conductivity	1:10 verd./ dil. in H ₂ O		5.9
Mindesthaltbarkeit minimum shelf life			15.06.2024
Lagerung (°C) storage			+15 bis/ to +30

Farbunterschiede von lot zu lot sind möglich.
The colour may change from batch to batch.

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 09.05.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.