

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS



IPG Chamber Cleaner	Kat.Nr./Cat.No. : 43399
IPG Chamber Cleaner	Lot/Contr.No. : P130057

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		gelbe, viskose Flüssigkeit yellow, viscous liquid	entspricht corresponds
Anwendung use		Reinigung cleaning procedures	
pH	1:10 verd./ dil. in H ₂ O	7.0 – 8.0	7.0
Leitfähigkeit (mS/ cm) conductivity	1:10 verd./ dil. in H ₂ O		5.9
Mindesthaltbarkeit minimum shelf life			15.06.2025
Lagerung (°C) storage			+15 bis/ to +30

Farbunterschiede von lot zu lot sind möglich.
The colour may change from batch to batch.

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 17.04.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.