

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS



SERVA PRiME™ Lightning Red	Kat.Nr./Cat.No. :	43402
	Lot/Contr.No. :	250184

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Komponenten components	SERVA PRiME™ Lightning Red DMSO, wasserfrei/ anhydrous Markierungspuffer/ Labeling Buffer Probenpuffer/ Loading Buffer		
Anwendungstest application test	SDS-PAGE	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			18.01.2026
Lagerung (°C) storage			+2 bis / to +8

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 27.02.2025

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.