

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

Visking® Dialysierschlauch	Kat.Nr./Cat.No. :	44114
Visking® dialysis tubing	Lot/Contr.No. :	130825

Parameter parameter	Spezifikation specification	Ergebnis result
Material	regenerierte Cellulose regenerated cellulose	entspricht corresponds
Flachbreite (mm) flat width	34	entspricht corresponds
Durchmesser (mm) diameter	21	entspricht corresponds
Ausschlußvolumen (Dalton) exclusion limit	12 000 – 14 000	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		07/ 2018
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 18.07.2013

Dipl.-Ing. (FH) Bernhard Göckel

Jakob Arndt

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.