

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

Servapor® Dialysierschlauch, MWCO 12 000 - 14 000	Kat.Nr./Cat.No. : 44148
Servapor® dialysis tubing, MWCO 12 000 - 14 000	Lot/Contr.No. : 240450

Parameter parameter	Spezifikation specification	Ergebnis result
Material material	regenerierte Cellulose regenerated cellulose	entspricht corresponds
Flachbreite (mm) Nominal dry flat width	80	entspricht corresponds
Durchmesser (mm) Nominal dry diameter	50	entspricht corresponds
Porendurchmesser (Å) pore diameter	ca. 25	entspricht corresponds
Ausschlußvolumen (Dalton) exclusion limit	12 000 – 14 000	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		04.03.2029
Lagerung (°C) storage		+15 bis/ to +30

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 22.05.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.