

**A N A L Y S E N Z E R T I F I K A T**  
CERTIFICATE OF ANALYSIS

<b>Serdolit® Rot p.A.</b>	<b>Kat.Nr./Cat.No. :</b>	<b>45540</b>
Serdolit® Red analytical grade	<b>Lot/Contr.No. :</b>	<b>121358</b>

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Aussehen</b> appearance		purpurvioletes Korn purple violet particles	entspricht corresponds
<b>Kapazität (eq/l)</b> capacity		min. 1.6	1.64
<b>Gegenion</b> counter-ion		H <sup>+</sup>	entspricht corresponds
<b>Korngröße (mesh)</b> particle size		16 – 50 (0.3 – 1.2mm)	entspricht corresponds
<b>Trocknungsverlust (%)</b> loss on drying	105°C; 3h	50.0 – 60.0	52.6
<b>Mindesthaltbarkeit</b> minimum shelf life			10/ 2017
<b>Lagerung (°C)</b> storage			+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date:** 11.12.2012

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**

This report has been computer generated and does not contain a signature.