

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Serdolit® Rot p.A.	Kat.Nr./Cat.No. :	45540
Serdolit® Red analytical grade	Lot/Contr.No. :	141106

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		purpurvioletes Korn purple violet particles	entspricht corresponds
Kapazität (eq/l) capacity		min. 1.6	1.70
Gegenion counter-ion		H ⁺	entspricht corresponds
Korngröße (mesh) particle size		16 – 50 (0.3 – 1.2mm)	entspricht corresponds
Trocknungsverlust (%) loss on drying	105°C; 3h	50.0 – 60.0	54.0
Mindesthaltbarkeit minimum shelf life			09/ 2019
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 21.10.2014

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.