

ANALYSENZERTIFIKAT
 CERTIFICATE OF ANALYSIS

Isopropanol für LC-MS	Kat.Nr./Cat.No. :	45636
Isopropanol for LC-MS	Lot/Contr.No. :	241063

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₃ H ₈ O	
Molekulargewicht molecular weight		60.10	
Aussehen appearance		klare, farblose Flüssigkeit clear, colourless liquid	entspricht corresponds
Gehalt (%) assay	GC	min. 99.95	100
Brechungsindex refractive index	20 °C	1.375 – 1.379	1.377
Azidität/ acidity (%)		≤ 0.0010	< 0.001
Alkalinität/ alkalinity (%)		≤ 0.0005	< 0.0005
Wasser/ water (ppm)	KF	≤ 300	127
Abdampfrückstand (ppm) residue on evaporation		≤ 2	< 2
Transmission (%)	220 nm	min. 64.0	68
	230 nm	min. 80.0	83
	260 nm	min. 98.5	98.5
HPLC Gradient (mAU)	254 nm	max. 2	< 1
Sensitive impurities (ppb) (reserpine)	Test LC-MS TIC (50 – 2000 m/z) ESI(+)	max. 100	< 100
Metalle/ metals (ppb)	Al, Fe, Na, Ca, Mg, K	max. 50 je/ each	< 50 je/ each
Mindesthaltbarkeit minimum shelf life			03.12.2027
Lagerung/storage (°C)			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.

Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 05.12.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.