

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

Colchicin Lösung sterilfiltriert	Kat.Nr./Cat.No. : 47252
Colchicin solution sterile filtered	Lot/Contr.No. : P240330

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		klare farblose Lösung clear colorless solution	entspricht corresponds
Konzentration (µg/ml) concentration	in PBS	10	10
pH		7.0 – 8.0	7.3
Osmolarität (mOsmol/l) osmolarity		270 – 300	280
Sterilität sterility		negativ negative	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			28.10.2026
Lagerung (°C) storage			-15 bis/ to -25

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 04.11.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.