

ANALYSENZERTIFIKAT
 CERTIFICATE OF ANALYSIS

Collagen R Lösung 0.2 %, steril	Kat.Nr./Cat.No. : 47254
Collagen R solution 0.2 %, sterile	Lot/Contr.No. : 240718

Parameter parameter	Spezifikation specification	Ergebnis result
Molekulargewicht molecular weight	ca. 300 000	
Aussehen appearance	viskose Lösung viscous solution	entspricht corresponds
Konzentration (%) concentration	2 mg/ ml in Essigsäure 0.1% 2 mg/ ml in acetic acid 0.1%	entspricht corresponds
pH	3.2 – 3.8	3.5
Sterilttest sterile test	kein Wachstum no growth	entspricht corresponds
Gelbildung gel formation	positiv positive	entspricht corresponds
Ursprung origin	Rattenschwanz rat tail	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		07.08.2027
Lagerung (°C) storage		+2 bis/ to +8

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control
Druckdatum/ printing date: 27.08.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.