

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

LB, Pulvermedium	Kat.Nr./Cat.No. : 48501
LB medium, powder	Lot/Contr.No. : P150300

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		elfenbeinfarbenes Pulver ivory-colored powder	entspricht corresponds
Löslichkeit solubility	20 g/l in H ₂ O	goldgelbe Lösung gold-yellowish solution	entspricht corresponds
pH	20 g/l in H ₂ O	7.0 ± 0.3	7.0
Osmolarität (mOsmol/l) osmolarity	20 g/l in H ₂ O		256
Zusammensetzung composition			
Tryptone (# 48647)		10.0 g/l	entspricht corresponds
Yeast extract (# 24540)		5.0 g/l	entspricht corresponds
NaCl (# 30183)		5.0 g/l	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			26.08.2018
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 28.08.2015

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.