

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

LB Agar, Pulvermedium	Kat.Nr./Cat.No. : 48502
LB Agar, powder	Lot/Contr.No. : P140301

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		hellelfenbeinfarbenes Pulver brightly ivory-colored powder	entspricht corresponds
Löslichkeit solubility	35 g/l in H ₂ O	goldgelbe Lösung gold-yellowish solution	entspricht corresponds
Zusammensetzung composition			
Tryptone (# 48647)		10.0 g/l	entspricht corresponds
Yeast extract (# 24540)		5.0 g/l	entspricht corresponds
NaCl (# 30183)		5.0 g/l	entspricht corresponds
Agar (# 11396)		15.0 g/l	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			07/ 2017
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control**Druckdatum/ printing date: 09.07.2014**

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.